




Concept
Medical

MagicTouch **SCB**


SIROLIMUS COATED BALLOON CATHETER



SIROLIMUS

MagicTouch **SCB**

SIROLIMUS COATED BALLOON CATHETER



Sirolimus Coated Polymer Free Balloon for Treatment of ISR, De-novo, Bifurcation and Small Vessel Disease Lesions.

- **CIRCUMFERNATIAL COATING**

Sirolimus Coating done under low pressure inflation for uniform coating and distribution of Drug

- **NANOLUTE TECHNOLOGY**

Encapsulation of sub micron sized Sirolimus particles in phospholipid carriers for effective transfer to deepest layer of vessel

- Minimal in-transit drug loss due to better re-wrap of balloon
- Presence of In-tissue Drug for treatment upto 120 days
- Programmed drug release ensure Sirolimus addresses therapeutic window
- Low tip entry profile for better trackability and deliverability
- Gold standard in SCB with 80000+ patients treated globally across CV disease spectrum
- **Highest studied SCB with IDE clearance**
A robust clinical trial program of more 6000+ patients studied globally across CV disease complexity